

# Rabbit Adrenal Medullary Cells

Cat. No. ARP0765,  $5 \times 10^5$  cells/vial

## Description

Research on the Rabbit Adrenal Medullary Cells is essential to the study of pheochromocytoma models, stress-induced catecholamine dysregulation, hypertension studies, autonomic nervous system dysfunction, and neuroendocrine tumor research. The adrenal glands are paired, triangular-shaped endocrine organs situated on the superior (upper) pole of each kidney. Each adrenal gland contains two main parts: the cortex and the medulla. The cortex produces steroid hormones, such as cortisol, aldosterone, and androgens, while the medulla secretes catecholamines, such as adrenaline and noradrenaline. These adrenal-gland-secreted hormones regulate essential physiological functions, including the aspects of metabolism, stress response, osmoregulation, and reproductive function. Dysfunction of the adrenal glands potentially results in conditions such as hypertension, metabolic imbalances, or endocrine disorders. The Rabbit Adrenal Medullary Cells are to be used with Rabbit Adrenal Medullary Cell Medium (Cat. No. ACM0765). This product is intended for laboratory in vitro use only. It is not intended for diagnostic, therapeutic, or clinical applications.

## Specification

Cell Type: N/A

Tissue/Organ: Adrenal gland

Disease: Normal

Species: *Oryctolagus cuniculus* (Rabbit)

Genetic Background: N/A

Markers: Neuron-Specific Enolase (NSE)

Symbols: RaAMC

## Shipping & Storage

Shipping condition: Frozen on dry ice.

Storage condition: Liquid nitrogen (LN<sub>2</sub>) cryopreservation.



## Intended Use

This product is intended for laboratory in vitro use only. It is not intended for diagnostic, therapeutic, or clinical applications.

## Culturing Guidance

Morphology: Fusiform, Irregular

Growth Mode: Adherent

Temperature: 37°C

Atmosphere: 5% CO<sub>2</sub>

## Unpacking and Storage Instructions

1. Visually inspect all packaging components for integrity and verify adequate dry ice.

If any damage is observed, notify Ascent Technical Support immediately.

2. Prioritize transfer to liquid nitrogen vapor phase storage system (-130°C or below).

Secondary option: -80°C mechanical freezer (short-term storage only).

Always maintain temperature strictly below -65°C.

## Disclaimer

Ascent Research endeavors to provide accurate and up-to-date product information. However, no warranties or representations are made regarding its completeness or reliability. References to scientific literature and patents are for informational purposes only, and the customer assumes sole responsibility for verifying their accuracy.

By accepting this product, the customer acknowledges and agrees to assume all risks associated with its receipt, handling, storage, disposal, and use, including compliance with all applicable safety and environmental regulations and precautions. Relevant laws, regulations, and ethical guidelines must be followed in conducting any research, modifications, or derivatives derived from this product.

This product is provided "AS IS", and except as expressly stated herein, Ascent Research disclaims all other warranties, express or implied. Under no circumstances shall Ascent Research, its affiliates, or representatives be liable for indirect, incidental, consequential, or punitive damages arising from the use of this material. While Ascent Research employs rigorous quality control measures, we shall not be held responsible for damages resulting from misidentification or misinterpretation of the provided materials.

# Copyrights

© 2025 Ascent Research. All rights reserved.

This document was last updated on June 20, 2025.

