

# Rabbit Fallopian Tube Smooth Muscle Cells

Cat. No. ARP0777,  $5 \times 10^5$  cells/vial

## Description

Research on the Rabbit Fallopian Tube Smooth Muscle Cells is essential to the study of tubal motility disorders, post-inflammatory tubal dysfunction, ovum transport impairment, and adenomyosis-related tubal spasm. The fallopian tubes are paired ducts that connect the ovaries to the uterus. Their primary function includes transporting oocytes, providing the site for fertilization, and delivering the zygote to the uterus for implantation. Obstruction or dysfunction of the fallopian tube can lead to infertility or an ectopic pregnancy. Additionally, the fallopian tubes are a common target for surgical sterilization procedures, such as tubal ligation, to prevent pregnancy. The Rabbit Fallopian Tube Smooth Muscle Cells are to be used with Rabbit Fallopian Tube Smooth Muscle Cell Medium (Cat. No. ACM0777). This product is intended for laboratory in vitro use only. It is not intended for diagnostic, therapeutic, or clinical applications.

## Specification

Cell Type: Muscle Cells

Tissue/Organ: Fallopian tube

Disease: Normal

Species: *Oryctolagus cuniculus* (Rabbit)

Genetic Background: N/A

Markers:  $\alpha$ -Smooth Muscle Actin ( $\alpha$ -SMA)

Symbols: RaFTSMC

## Shipping & Storage

Shipping condition: Frozen on dry ice.

Storage condition: Liquid nitrogen (LN<sub>2</sub>) cryopreservation.

## Intended Use



This product is intended for laboratory in vitro use only. It is not intended for diagnostic, therapeutic, or clinical applications.

## Culturing Guidance

Morphology: Fusiform, Irregular

Growth Mode: Adherent

Temperature: 37°C

Atmosphere: 5% CO<sub>2</sub>

## Unpacking and Storage Instructions

1. Visually inspect all packaging components for integrity and verify adequate dry ice.  
If any damage is observed, notify Ascent Technical Support immediately.
2. Prioritize transfer to liquid nitrogen vapor phase storage system (-130°C or below).  
Secondary option: -80°C mechanical freezer (short-term storage only).  
Always maintain temperature strictly below -65°C.

## Disclaimer

Ascent Research endeavors to provide accurate and up-to-date product information. However, no warranties or representations are made regarding its completeness or reliability. References to scientific literature and patents are for informational purposes only, and the customer assumes sole responsibility for verifying their accuracy.

By accepting this product, the customer acknowledges and agrees to assume all risks associated with its receipt, handling, storage, disposal, and use, including compliance with all applicable safety and environmental regulations and precautions. Relevant laws, regulations, and ethical guidelines must be followed in conducting any research, modifications, or derivatives derived from this product.

This product is provided "AS IS", and except as expressly stated herein, Ascent Research disclaims all other warranties, express or implied. Under no circumstances shall Ascent Research, its affiliates, or representatives be liable for indirect, incidental, consequential, or punitive damages arising from the use of this material. While Ascent Research employs rigorous quality control measures, we shall not be held responsible for damages resulting from misidentification or misinterpretation of the provided materials.

## Copyrights

© 2025 Ascent Research. All rights reserved.

This document was last updated on June 20, 2025.

